
Cystocele Repair Patch®

3 mesi dall'implanto - 3 months post implant - 3 meses después del implante

BioRipar®
Sede di impianto
Site of implant
Sitio del implante

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DESCRIPTION

Pelvic organic support, described by J. DeLancey, (J.O.L. DELANCEY 1994) identifies in the level II the anchorage of the medial portion of the vaginal walls through the pubo-cervical fascia, which is fixed to the tendinous arch of the endopelvic fascia. Therefore, the anatomic integrity is the main mandatory condition to the functional restore of pelvic organs.

STANDARD PELVIC SURGERY

Standard pelvic floor surgery, main phases:

- I.** Anterior vaginal wall incision «reverse T» form the vesical neck to the portio, bilaterally, or to the vaginal vault, in case of previous hysterectomy
- II.** Pubic-vesical fascia separation or "Alban's fascia" from the vaginal wall to the ischio-pubic branch
- III.** Pubic-vesical fascia duplication through monofilament long term absorbable suture Cal. 2/0 (Assuofil Monofilament FU395M), separate stitch technique
- IV.** Vaginal wall enclosure
- V.** Vaginal Breach closure through continuous stitch technique with barbed long term absorbable suture (Filbloc FQ07GHHAD).

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GUIDED REPAIR PATCH

The integration in the host tissue three months after implantation, shows an advanced integration, with evident deposit of neo-collagen and angiogenesis. The process will finish with a physiological regeneration of a new functional structure with increased tensile strength, useful to repair the anatomic deficit of cystocele, the compartment with the highest relapse rate in pelvic surgery.

FUNCTIONAL ANATOMIC REPAIR

CRPatch Cystocele Repair Patch has a specific design developed to optimize the surgical strategies oriented to the functional anatomy in terms of resolution of the specific defect. The shape has been designed to recreate the anatomy of DeLancey's level II, reinforcing specifically the damaged tissues, fixing the patch's wings to the tendinous arch of the endo-pelvic fascia. The sub-urethral positioning and the CRPatch's medial body perform the repair of the central defect of the cystocele, the lateral wings are used to solve the lateral defect by functional anchorage. The specific procedure is the same compared to the standard surgery for the first 3 steps, then, the paravesical spaces are prepared, for a depth of about 3 cm, towards the tendinous arch, allocating the lateral wings in this site, immediately after the sub-urethral positioning. The fixation to the pubo-cervical fascia is performed with Cal. 2/0 monofilament long term absorbable suture (Assufil Monofilament FU395M) with separate stitch technique.

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