



ASSUT EUROPE
SUTURE CHIRURGICHE



More than surgery




COVID-19

ANTIGEN RAPID TEST
NASOPHARYNGEAL SWAB

SIMPLE

EASY

ACCURATE



“SARS-CoV-2 Antigen Detection Kit” is an in vitro immunochromatographic assay for qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in nasopharyngeal (NP) swab specimens.

It is intended to aid in the diagnosis of SARS-CoV-2 infections. It does not differentiate between SARS-CoV and SARS-CoV-2.


SPECIMEN COLLECTION



Use the nasopharyngeal swab supplied in the kit

1. Carefully insert the swab into the nostril of the patient, reaching the surface of posterior nasopharynx
2. Swab over the surface of the posterior nasopharynx. Rotate the swab several times
3. Withdraw the swab from the nasal cavity

SAMPLE PREPARATION PROCEDURE



1. Insert the test extraction tube into the workstation in this product. Make sure that the tube is standing firm and reaches the bottom of workstation


2. Add 0.3mL (about 10 drops) of the sample extraction buffer into the extraction tube

3. Insert the swab into the extraction tube which contains 0.3mL of the buffer extraction


4. Roll the swab at least 6 times while pressing the head against the bottom and side of the extraction tube

5. Leave the swab in the extraction tube for 1 minute

6. Squeeze the tube several times with finger from outside of the tube to immerse the swab. Remove the swab. The extracted solution will be used as test sample.



INTERPRETATIONS OF RESULTS



Positive (+): Both of T and C lines are appeared in 15minutes after the sample added.

Negative (-): C line is appeared while no T line appeared in 15 minutes after the sample added.

Invalid: FAs long as the C line does not appear, it indicates that the test result is invalid, and should retest with another test card.

Wait for the colored band (s) to appear.

Read the result in 15 minutes. Do not interpret the result after 20 minutes.

KIT COMPONENTS

- Test cassette disposable
- Extraction tube and dropper tip
- Nasopharyngeal swab disposable
- Buffer

PRODUCT CODE	DESCRIPTION	BOX
COVID-19/AGA/40	SARS-CoV-2 Antigen Detection Kit	BOX of 40



More than surgery



ITALIA

Stabilimento di Produzione

Zona Industriale
67062 Magliano dei Marsi (AQ), Italy
Tel. +39 0863 517956/515000
Fax +39 0863 570084/515209

Ufficio Legale

Via G. Gregoraci, n.12
00173 Roma
Tel. +39 06 72677348
Fax +39 06 72675380

FRANCE

86-114 Avenue Louis Roche
Batiment D Hall-3 92230
Gennevilliers
Tel. +33 147985294
Fax +33 147985892

ESPAÑA

Poligono Landaben
Calle C, nave 1y2 - 31012
Pamplona (Navarra)
Tel. +34 948122274
Fax. +34 948120842

GERMANY

Prager Strasse 34
04317 Leipzig

UK

Unit 12a Barncliffe Mills
Shelley HD8 8LU Huddersfield
Tel. +44 1484602222

USA

55 NE 5th Avenue, Suite 501
33432 - 4093 Boca Raton FL
United States

RUSSIAN FEDERATION

c/o Studio Marzona
Rozhdestvenskiy B-r 10 Bld. 7
107031 Moscow
Tel/Fax +7 962 9862804

BRASIL

Rua Professor Alfredo Gomes, 14
CEP 22551-080
Botafogo Rio De Janeiro (RJ)
Tel.: 0055 21 22661943
Fax: 0055 21 22663539

ALGERIE

Zone D'Activité Dar El Beida
Lot n° 66, Bt 12
Alger

TUNISIE

Rue 7316 Villa n°5 - El Menzah 9B
1013 Tunis, Tunisie
Tel.: 00216 71875740
Fax: 00216 71880541

www.assuteurope.com

info@assuteurope.com



©2020 Assut Europe S.p.A. All rights reserved
CE 0373 - UNI EN ISO 9001:2015
UNI EN ISO 13485:2016
UNI EN ISO 14001:2015